trophon[®] EPR

User Manual

trophon[®] EPR

User Manual

Read this manual before operating the trophon® EPR to determine the correct procedures.

For further information, contact your customer service representative or visit the Nanosonics website.

All technical specifications and system approvals are listed in Appendix 1.

©2015 Nanosonics Limited. All rights reserved.

The content of this manual is correct at the time of product purchase.

trophon®, NanoNebulant® and Sonex-HL® are registered trademarks of Nanosonics Limited.

NanoNebulant is the product name of the trophon disinfectant used in all regions where trophon is available for sale, with the exception of US and Canada.

Sonex-HL is the product name of the trophon disinfectant in the US and Canada.

Nanosonics' technologies are protected by a combination of patents, trademarks, and exclusive proprietary rights <http://www.nanosonics.com.au/nanosonics-patent>.

Your sales representative is:

Attach Business Card or information sticker/stamp here.

 0197 Conforming to MDD 93/42/EEC; certified by Notified Body TUV Rheinland

Conforming to RoHS Directive 2011/65/EC

 Conforming to WEEE Directive 2002/96/EC under Article 7 Recovery

www.nanosonics.eu • www.nanosonics.com.au • www.nanosonics.us


Table of Contents


PART A – WARNINGS, INSTRUCTIONS AND INTRODUCTION	4
SECTION A1: IMPORTANT WARNINGS, LABELS and SYMBOLS	4
A1.1 Labels and Symbols	4
A1.2 Instructions	5
A1.3 Warnings	5
SECTION A2: INTRODUCTION TO THE TROPHON EPR	6
A2.1 Intended Use	6
A2.2 Disinfection Process	6
A2.3 Disinfection History Log	6
A2.4 Validated Probes and Disinfectants	6
A2.5 Training	6
A2.6 Environment and User Profile	6
A2.7 Accessories	6
PART B – OVERVIEW OF TROPHON FEATURES AND INSTALLATION GUIDE	8
SECTION B1: TROPHON FEATURES	8
SECTION B2: INSTALLATION GUIDE	9
B2.1 Positioning your trophon	9
B2.2 Cable Management System Installation	9
B2.3 Powering On	10
B2.4 Basic Settings	10
B2.5 Warm-up Cycle	10
PART C – ROUTINE USE, MAINTENANCE AND CARE	11
SECTION C1: LOADING THE DISINFECTANT CARTRIDGE	11
SECTION C2: ROUTINE HIGH LEVEL DISINFECTION CYCLE	11
C2.1 Preparing the Probe	11
C2.2 Inserting the Chemical Indicator	11
C2.3 Positioning the Probe	11
C2.4 Closing the Chamber Door	13
C2.5 Disinfecting the Probe	13
C2.6 Removing the Probe	14
C2.7 Sleep Mode and Shutdown Procedures	14
SECTION C3: PURGE CYCLE	15
C3.1 When to Run a Purge Cycle	15
C3.2 How to Manually Run a Purge Cycle	15
C3.3 Transporting the trophon	15
SECTION C4: INCOMPLETE OR FAILED CYCLES	15
C4.1 Mains Power Failure	15
C4.2 Cycle Fault	15
C4.3 Manual Door Lock Override	16
SECTION C5: ROUTINE CARE AND MAINTENANCE	16
C5.1 Regular Cleaning	16
C5.2 Service	16
SECTION C6: DISPOSAL OF TROPHON	16
PART D – TROUBLESHOOTING	17
PART E – SERVICE AND WARRANTY	18
APPENDIX 1: TROPHON EPR TECHNICAL SPECIFICATIONS, N00010, N00020, N03000 & N04000	18
APPENDIX 2: PRODUCT WARRANTY TERMS AND CONDITIONS	19

PART A – WARNINGS, INSTRUCTIONS AND INTRODUCTION

SECTION A1: IMPORTANT WARNINGS, LABELS and SYMBOLS

A1.1 Labels and Symbols

	Caution		Corrosive
	Consult instructions for use		Warning
	Start (of action)		Single Use Only
	Fragile / Handle With Care		Warning: Hot Surface
	Do not disassemble		Dangerous Voltage
	Separate collection for electrical and electronic equipment		Keep Dry
	Keep Out of Direct Sunlight		Expires (year and month)
	Batch Number		Product Number
	Cannot be transported by air freight		This Way Up
	UN 2014		Wear Gloves when handling cartridges
	Oxidizer – 5.1		Corrosive – 8
	Legal Manufacturer		Date of Manufacture
	trophon EPR Storage and Transport Conditions. Temperature range: -20°C to +60°C / -4°F to +140°F		Environmental Specifications. trophon EPR operating temperature range: 17°C to 27°C / 62.6°F to 80.6°F

Probe Positioning Guide

A1.2 Instructions

Read these instructions before using the trophon EPR:

- Safety Data Sheet (SDS) enclosed with the trophon.
- Occupational or Workplace Health and Safety Guidelines (OH&SG, OSHA, WHS) for your institution for lifting, spills etc.
- Chemical Indicator Instructions For Use (IFU)
- trophon disinfectant cartridge IFU
- Instructions for trophon EPR accessories.
- Probe manufacturer's instructions.

If you do not follow instructions:

- Burns, bleaching, electric shock or other injury may occur.
- High level disinfection may not be achieved.
- Residual disinfectant may remain on the probe, which may cause injury when removing the probe.
- Equipment damage may occur.

A1.3 Warnings

Hot Temperatures

- Do NOT touch surfaces in the internal chamber. These are hot and can cause burns.
- To avoid probe damage, ensure the probe is correctly positioned in the chamber. (See section C2.3 for correct positioning of probe.)

Malfunctions

- Do NOT attempt to open the chamber door during a cycle, power failure or system malfunction.
- All repairs must be carried out by trained personnel.

Transporting the trophon

- The trophon weighs approximately 17kg (38lb).
- The packaged trophon weighs approximately 21kg (46lb)
- If your trophon has been used, purge the trophon before moving, to remove hydrogen peroxide. (See section C3).

Electrical Device

- Use the power cable supplied with the trophon. Connect to an earthed power outlet with the correct voltage and frequency as on the product and specified in Appendix 1. Incorrect voltage can cause damage.
- Spilled fluid can result in electrical shock. Avoid spilling fluids onto or around the trophon. Do not immerse any parts in liquid.
- Do not attempt to access the internal mechanics. This may result in electrical shock.

Protective Wear and Spills

- Wear clean disposable gloves when running the trophon EPR and handling:

- disinfectant cartridges, as temporary bleaching and/or irritation of the skin may occur if no gloves are worn
- probes before and after a disinfection cycle
- emptying the waste container.
- Wear appropriate personal protective equipment (PPE) when managing spills.
- Never return spills to original cartridges for re-use.

SECTION A2: INTRODUCTION TO THE TROPHON EPR

A2.1 Intended Use

The trophon EPR together with the disinfectant cartridge is a high level instrument-grade disinfectant system designed to high level disinfect validated ultrasound probes according to the specified processes outlined in this manual. Do **NOT** use this device for any application other than its expressed purpose.

The trophon EPR is **NOT** intended to reprocess single-use devices or pre-clean ultrasound probes.

Chemical Indicator use is required with every disinfection cycle.

A2.2 Disinfection Process

At the start of the high level disinfection cycle, the trophon EPR creates an aerosol of concentrated hydrogen peroxide. This is distributed over the exposed surface of the probe providing thorough, high level disinfection of the shaft and the handle of the probe. The trophon then breaks down the hydrogen peroxide into small amounts of oxygen and water which is evaporated into the atmosphere. During a purge, oxygen is vented into the atmosphere and water is collected in the waste container inside the trophon (maximum capacity 150ml).

The contact conditions are fixed cycle parameters that the end user cannot modify.

A2.3 Disinfection History Log

A disinfection history log can be downloaded from your trophon EPR and printed onto to a PC or laptop via trophon Connect software (optional accessory - available for certain trophon variants only) or provided by your service centre during scheduled preventative maintenance services.

A2.4 Validated Probes and Disinfectants

For details of probes that can be used in the trophon EPR refer to:

- trophon EPR Validated Probes List at the Nanosonics website.

Use only trophon disinfectant cartridges when high level disinfecting with the trophon EPR.

A2.5 Training

Before setting up or using your trophon EPR, ensure that all users:

- Are trained in safety procedures and potential hazards, as outlined in this manual.
- Visit **Trophon-EPR/Login** on the Nanosonics website to complete the online training module and receive trophon EPR certification.

A2.6 Environment and User Profile

The trophon EPR is designed for use in healthcare facilities to high level disinfect ultrasound transducers under the control of trained healthcare professionals.

The disinfection cartridge, chemical indicator, and trophon EPR system are designed to be used with minimal personal protective equipment (gloves only), and in standard workplace or clinical settings, including at the patient point-of-care. Special ventilation and other safety precautions are not required when used as per the manufacturer's instructions.

A2.7 Accessories

Accessories* designed for use with the trophon EPR include:

- trophon Connect
- trophon Printer & Label Roll

- trophon Wall Mount
- trophon Cart
- trophon Logbook
- trophon Curved Probe Positioner (CPP)
- Clean Ultrasound Probe Cover
- trophon Printer Wall Mount
- trophon Printer Cart Mount

Only use trophon approved accessories or the trophon may be ineffective.

***Some accessories are available for certain trophon variants only. Contact your trophon representative for more information.**

PART B – OVERVIEW OF TROPHON FEATURES AND INSTALLATION GUIDE

SECTION B1: TROPHON FEATURES

Front

Figure 1

Back


Figure 2

Chamber


Figure 3

Left Side


Figure 4

Right Side


Figure 5

1. User display screen
2. Start button
3. Soft key buttons
4. Chamber door handle
5. Adjustable feet (for levelling)
6. Chamber door (opened)
7. Probe in correct position
8. Cable clamp
9. Gland seal
10. Power switch
11. Power socket
12. Serial port
13. Cartridge door- **Do NOT force open.**
Cartridge door will automatically open when required.
14. Waste container

SECTION B2: INSTALLATION GUIDE

B2.1 Positioning your trophon

1. Ensure the surface can support the weight and allows adequate airflow (see Appendix 1).
2. Ensure the area around your trophon is free from other equipment and clutter. Position as shown in Figure 6 to ensure access to all features.

3. Turn the back feet either clockwise or anti-clockwise to level (see Figure 2). Check feet are not loose or missing.

The trophon EPR should be placed at a height from the floor level to accommodate a range of user heights. Refer to Figure 7 for a guide to an ergonomically safe work zone for most users.


Figure 7

B2.2 Cable Management System Installation

The Cable Management System (CMS) holds the ultrasound cable away from the chamber door, and is storage for the cable during high level disinfection.

To install the CMS:

1. Open the chamber door.
2. Peel back the protective strip from the double sided tape on the bottom of the CMS (see Figure 8).
3. Align and slide the two hooks on the front of the CMS into the slots in the trophon chamber (see Figure 9-1) until it cannot slide any further.
4. Push the CMS down on top of the trophon to secure it in place (see Figure 9-2).


Figure 8

Figure 9

B2.3 Powering On

1. Attach the power cable supplied to the trophon power socket.
2. Connect to mains power and switch on.
3. Turn on the power switch, located on the side of the trophon.

NOTE: To maximise the life of the disinfectant cartridge, keep your trophon switched on at all times, except when the machine needs to be moved.

B2.4 Basic Settings

For all settings:

1. Press the soft key underneath *Menu* on the screen.
2. Using the soft keys under the LCD:
 - Scroll to *Setup* and press *OK*.
 - Press *OK* after each setting has been made.

Date and Time

- Scroll to *Set Date and Time* and follow the onscreen messages.

Language

- Scroll to *Language* and follow the onscreen messages.

Alarm Settings

- Scroll to *Alarm Settings* and press *OK*.
- Change the alarm repeat or volume settings by following the onscreen messages.
- If alarm repeat is enabled, the alarm will repeat every 30 seconds for the following situations:
 - a. The probe is left in the trophon after a HLD cycle.
 - b. An error message is displayed and requires a response.

Contrast

- Scroll to *Contrast* and follow the onscreen messages.

B2.5 Warm-up Cycle

1. The warm up cycle prepares the trophon for operation and will begin automatically when the machine is powered on or restarted from sleep.

Screen Message	Approximate Warm Up Time (minutes)
<i>Quick Warm Up</i>	< 2
<i>Warming Up</i>	2–30
<i>Extended Warm Up</i>	> 30

2. Instructions may appear on the screen during this time, which should be followed, such as:
 - *Close Chamber Door*
 - *Cartridge empty. Replace Cartridge Now?*
 - *Remove Probe from Chamber*
3. When warm up is complete, the screen message will read *Load Probe and Indicator*.
4. The machine is now ready for use.

PART C – ROUTINE USE, MAINTENANCE AND CARE

SECTION C1: LOADING THE DISINFECTANT CARTRIDGE

A trophon disinfectant cartridge needs to be inserted into the trophon EPR before a high level disinfection cycle can commence. Refer to the trophon disinfectant IFU for detailed instructions on how to insert and/ or replace trophon disinfectant cartridges.

SECTION C2: ROUTINE HIGH LEVEL DISINFECTION CYCLE

C2.1 Preparing the Probe

Wearing gloves, clean and dry the probe BEFORE commencing the high level disinfection process as per the probe manufacturer's instructions.

C2.2 Inserting the Chemical Indicator

A trophon Chemical Indicator must be used for each disinfection cycle. No other chemical indicators are approved for use in the trophon EPR. Each Chemical Indicator may only be used once. Place the trophon Chemical Indicator into the holder. Refer to the trophon Chemical Indicator IFU. You may then position the probe in the chamber.

C2.3 Positioning the Probe

1. When the trophon is ready, screen message *Load Probe and Indicator* will be displayed.
2. Open chamber door.
3. Two clamps hold the probe securely in the chamber. See Figure 10.

The probe has a short sleeve at the back of the handle, covering the electrical cable. This is the Probe Gland. See Figure 10(a).

Wearing gloves insert the probe correctly in the trophon and ensure the probe gland is held by the soft clamp. See Figures 10–12.

Ensure the probe is correctly positioned in the chamber. The probe must not contact the chamber wall and must be positioned above the embossed line at the bottom of the chamber.


Figure 10


Figure 11


Figure 12

4. If a cable management system has been installed on your trophon (see section B2.2), secure the external portion of the cable and connector by:

- Placing the connector carefully inside the tray.
- Coiling the cable neatly and safely around the cable holder on the side of the tray assembly.

NOTE: Incorrect positioning of the probe may result in:

- High-level disinfection will not be achieved during the trophon EPR disinfection cycle
- Residual disinfectant remaining on the probe surface. This may lead to temporary bleaching and/ or irritation of the skin if no gloves are worn.
- Damage to the probe may occur if in contact with the chamber wall

NOTE: Curved probes must be correctly inserted in the trophon. See Figure 13. The curved probe positioner (CPP) accessory is required to assist with this positioning for certain probes. Refer to the trophon EPR Validated Probes List.


Figure 13

C2.4 Closing the Chamber Door

- The door will automatically lock at the start of a high level disinfection cycle.
- If the door is not properly closed, screen message: *Close Chamber Door* will be displayed.

C2.5 Disinfecting the Probe

The next screen message: *Is the Probe Clean and Dry?*

1. If the probe has been pre-cleaned and dried, press *Yes* using the soft key button. If *Yes* is selected, you will be prompted by the screen message: *Press start to begin*.
2. If the probe has NOT been pre-cleaned and dried, press *No*. Remove the probe, clean and dry the probe and position again as directed in section C2.3.
3. Once the probe is pre-cleaned and dried, select *Yes* and press the **START** button to initiate the cycle or *Cancel* to unlock the machine and remove the probe.
4. The progress of the disinfection cycle is indicated on the LCD screen message: *Disinfecting*.

5. The High Level Disinfection cycle will take 7 minutes to complete.

WARNING: In the unlikely event mist is visibly escaping from the chamber, remain at a distance from the trophon until completion of the operating cycle and the mist is no longer visible. Do not come into direct contact with the mist. Contact your customer service representative (see Part D – Troubleshooting).

C2.6 Removing the Probe

1. When the cycle has been completed, the trophon will sound an audible alarm.

NOTE: After completion of a successful high level disinfection cycle, the ultrasound probe and chamber may have surface temperatures up to 45°C/ 113°F and 60°C/ 140°F respectively. Care should be taken not to touch the chamber. The probe will feel warm to touch and is safe for handling and use with gloves on.

2. Follow the on screen instructions:

Screen message 1: Cycle Complete. Wear Gloves. Remove and Wipe Probe.

Screen message 2: Attention: Wear Gloves and Wipe Probe. (This message indicates that some hydrogen peroxide may not have been broken down into its constituents, oxygen and water and extra care should be taken when removing the probe).

3. Wearing clean gloves, immediately remove the used Chemical Indicator from the trophon and verify the colour change against the chart on the Chemical Indicator carton. Record the result using the printer or logbook accessories if required.
4. Discard used Chemical Indicator.
5. Remove the probe carefully using minimal contact after the cycle is complete. Avoid touching the probe against the chamber's hot surface. Avoid compromising the probe's disinfected surface before use.

NOTE: After HLD cycle completion, the trophon performs a rapid cooling cycle until the probe is removed from the chamber to prevent overheating of the probe. If the probe is not removed immediately this will increase the warm-up time required by the subsequent cycle. It is therefore recommended to remove the probe as soon as possible after HLD cycle has been completed.

6. Wipe the probe with a clean, low lint, absorbent, single-use, dry cloth/ wipe. Visually inspect the probe and ensure any disinfectant residue present is removed.
7. If a pass was verified by the Chemical Indicator colour AND the trophon screen displayed *cycle complete*, the HLD has been successful – proceed to step 8 below. If one or both of these items did not occur, repeat steps from section C2.1.
8. Discard gloves.
9. Close the chamber door.
10. The probe is now ready for use or storage.

C2.7 Sleep Mode and Shutdown Procedures

- If the trophon is not used for 120 minutes or a probe has been left inside the trophon for an extended amount of time, it will automatically enter sleep mode in order to save power. To restart the trophon from sleep mode press *Restart*.
- While the trophon is in sleep mode it will perform self-maintenance functions and will display the messages: *Warming Up* or *System Refresh*. Please do not switch the trophon off during these processes.
- System refresh during sleep will only occur for low use customers and does not impact the number of disinfection cycles that can be performed per cartridge. This process will typically take 13 minutes.

SECTION C3: PURGE CYCLE

The purge cycle removes any remaining disinfectant from the cartridge as well as inside the trophon and converts the disinfectant into oxygen and water.

C3.1 When to Run a Purge Cycle

- When an onscreen message states that your trophon requires purging. This will occur upon disinfectant expiry at 30 days after insertion. Follow the onscreen instructions to purge.

Manually initiate a purge cycle:

- When the trophon detects an error that requires a service call.
- Before lifting or moving the trophon.
- For troubleshooting purposes when directed by service support only

C3.2 How to Manually Run a Purge Cycle

NOTE: Once the purge cycle has started it may be paused but it cannot be cancelled. Do not switch the trophon off during purging as this will restart the purge cycle. Do NOT attempt to open the cartridge door during the purge cycle.

1. To manually initiate a purge cycle:
 - Wearing gloves, ensure the empty waste container is fully inserted into the trophon.
 - Select *Menu*, scroll to *Purge* and select *OK*.
 - Confirmation screen message: *Remove All Disinfectant from System?*
 - A purge cycle will commence automatically once *Yes* is selected.
2. The purge cycle will typically take 35 minutes (maximum of 100 minutes).
3. When purging is complete, put on gloves and follow the screen messages: *Remove and Empty Waste Container* and *Load Waste Container*.
4. The next screen message to follow: *Cartridge Empty. Replace Cartridge Now?*
5. If *OK* is selected, the cartridge door will automatically open after *Safety Check* is performed and the empty cartridge can be removed. Do not force the cartridge door open.
6. Follow the next screen message: *Replace Cartridge and Close Cartridge Door*.

C3.3 Transporting the trophon

- Before transporting the trophon EPR, you must purge the disinfectant, switch off the trophon at the power switch and unplug the trophon from the mains
- Do not move excessively or drag the trophon across the mounting surface.
- Keep trophon upright at ALL times even during transportation.

SECTION C4: INCOMPLETE OR FAILED CYCLES

This section describes the most common situations in which a cycle has not been completed satisfactorily and the required actions to take (see also Part D – Troubleshooting).

C4.1 Mains Power Failure

If the mains power supply to the trophon is lost while in operation, the current cycle will not complete.

- Once power is restored follow the onscreen messages to remove the probe safely from the trophon.
- Discard the used Chemical Indicator and replace with a new one. Repeat the disinfection cycle.
- If power cannot be restored and the probe is urgently required, follow section C4.3.

C4.2 Cycle Fault

If a problem occurs during or at the end of the cycle, a cycle fault will be detected. Follow the onscreen messages to fix the fault and repeat the disinfection cycle.

In case of a repeated fault or serious malfunction, contact your customer service representative citing the error message shown on the screen. Do NOT attempt to use the trophon or the probe.

C4.3 Manual Door Lock Override


Use **ONLY** when the probe is locked in the chamber and must be **urgently** retrieved for use.

WARNING: There may still be disinfectant in the chamber and chamber surfaces may still be hot. Gloves must be worn to avoid contact with disinfectant.

By turning off the trophon and turning it back on, the trophon will attempt to recover. At this time if the probe still cannot be removed then the following may be undertaken:

1. Turn off the trophon.
2. Insert the tip of a screwdriver into the slot behind the chamber door handle gently until it stops. See Figures 14 and 15.
3. Lift the screw driver in an upwards motion.
4. The door will unlock and the probe can be removed.

NOTE: The probe is **NOT DISINFECTED** and **CANNOT** be used until it has completed a successful disinfection cycle or been disinfected by an alternative method.

SECTION C5: ROUTINE CARE AND MAINTENANCE

C5.1 Regular Cleaning

1. Do NOT submerge the trophon, or pour liquids over it.
2. Keep the trophon level and upright at all times.
3. Keep the power socket completely dry (see Figure 4).
4. For cleaning, when cool, wipe inside the chamber and the outside of the trophon with a cloth moistened with a mild, soapy solution, until all surfaces are visibly clean.
5. For disinfecting, wipe all accessible surfaces of the trophon with an Isopropanol or Quat Wipe.

C5.2 Service

Your trophon requires annual preventative maintenance as well as a major service at 5,000 cycles.

Once the service interval of 12 months or 5000 cycles is displayed on the screen, please contact your customer service representative to arrange a service.

SECTION C6: DISPOSAL OF TROPHON

Contact your distributor or Nanosonics. They will advise you of the collection point for recycling electrical and electronic equipment.

PART D – TROUBLESHOOTING

Use this table to diagnose and resolve basic problems. If a probe is present and you need to retrieve it urgently, follow section C4.3.

If the problem persists, contact your trophon representative.

Symptom	Check for the following:
1. There is no power to the trophon. 2. The screen is blank.	The trophon is completely plugged in and switched ON at both the trophon power switch and at the wall outlet. The correct power cable for your region is used.
3. The chamber door will not open.	There is power to the trophon. There is no disinfection, warm up or purge cycle running.
4. The chamber door will not close.	The probe is validated for use in the trophon EPR – see section A2.1. The probe is loaded correctly. The chamber door is not in the “locked position”: if in the ‘locked’ position, see Symptom 5 below.
5. The chamber door is open and locked.	If the trophon is in sleep mode, select <i>Restart</i> on the LCD; Otherwise, power the trophon OFF and back ON using the power switch (see Figure 4).The trophon should then automatically unlock the chamber door. If the above are not successful, power OFF the trophon and follow section C4.3. Close the chamber door before powering the trophon EPR ON again.
6. The cartridge door does not open.	There is power to the trophon. All cycles are complete. The waste container is empty and fully inserted. There is nothing obstructing the cartridge door from opening. Note, the cartridge door is automatic, and should not be forced open.
7. The cartridge door will not close.	Correct cartridge type has been inserted. Cartridge is correctly positioned. Cartridge lid has been removed.
8. The probe will not sit correctly in the chamber.	The probe is validated for use in the trophon EPR – see section A2.1. The probe is loaded correctly. Check whether a CPP accessory is required, on the Validated probes list.
9. The cycle will not start.	The probe is loaded correctly. The chamber door is closed. Confirm probe is clean & dry before pressing start button.
10. The trophon is beeping.	Probe inside the trophon which needs to be removed. The chamber door is open. On screen error message.
11. Liquid is leaking from the trophon.	 WARNING: Any fluid leaking from the trophon may contain hydrogen peroxide. If liquid or mist is seen coming from the trophon at any time: <ul style="list-style-type: none"> • Do not come into contact with the mist or liquid. • Wear appropriate PPE. • Ensure area is well ventilated. • Allow the trophon to complete the cycle. • Turn off the trophon and remove the power cord. • Contact your customer service representative.
12. The trophon is failing numerous cycles.	Record any error codes and colour of the chemical indicator and contact your customer service representative.
13. The trophon takes too long to warm up between cycles	Ensure that the probe is removed immediately at the end of each cycle.

PART E – SERVICE AND WARRANTY

Contact your Nanosonics representative if you have any questions about:

- The trophon EPR or accessories
- The warranty

Each trophon EPR has a comprehensive warranty against defects in material and workmanship for 12 months from the date of delivery.

Service Schedule: When due, this will be indicated prior to commencing a disinfection cycle by the screen message: *Device Service Due*, at weekly intervals until a service is performed. Service intervals can also be accessed via the *System Information*, located in the trophon EPR LCD menu.

NOTE: Use only authorised service personnel to service the trophon EPR.

Modifying the trophon without authorisation will void your warranty.

APPENDIX 1: TROPHON EPR TECHNICAL SPECIFICATIONS, N00010, N00020, N03000 & N04000

N00010 Electrical Specification	Rated mains input voltage: 120V AC Rated mains input current: 5Amp, 50/60Hz Mains Inlet: IEC type C13 Equipment must be connected to an earthed outlet using power cable supplied with the trophon.
N00020 Electrical Specification	Rated mains input voltage: 230V AC Rated mains input current: 5Amp, 50/60Hz Mains Inlet: IEC type C13 Equipment must be connected to an earthed outlet using power cable supplied with the trophon.
N03000 Electrical Specification	Rated mains input voltage: 230V AC Rated mains input current: 5Amp, 50/60Hz Mains Inlet: IEC type C13 Equipment must be connected to an earthed outlet using power cable supplied with the trophon.
N04000 Electrical Specification	Rated mains input voltage: 100V AC Rated mains input current: 6A, 50/60Hz Mains Inlet: IEC type C13 Equipment must be connected to an earthed outlet using power cable supplied with the trophon
Data port	Type: RS232 Connector: DB9 socket
Environmental Specification	Operating temperature range: 17 to 27°C / 62.6°F to 80.6°F Enclosure rated to IP20
Storage and Transport Conditions	Temperature range: -20°C to +60°C / -4°F to +140°F
Physical Characteristics	Weight of unpacked trophon: 17kg Dimensions of trophon: 485mm high x 345mm wide x 320mm depth
Electromagnetic Compliance	The trophon has been tested and found to comply with the limits for emission (electromagnetic Interference) pursuant to EN61326-1:2013 (CISPR 11 Group 1 Class B limits)

APPENDIX 2: PRODUCT WARRANTY TERMS AND CONDITIONS

Terms

This warranty is given by Nanosonics Limited ABN 11 095 076 896 of 14 Mars Road, Lane Cove, NSW 2066 Australia (**Nanosonics**).

Nanosonics warrants to the customer that the trophon EPR is free from defects in material and workmanship that materially affect its function under normal use and service for a period of 12 months commencing upon the date of purchase (**warranty period**).

You may have statutory rights in relation to the trophon EPR and these are not affected by this warranty.

Exclusions

This warranty does not apply in the following circumstances (regardless of how those circumstances arise):

- a. where the trophon EPR has not been used, handled, installed, stored, cleaned and serviced in accordance with the relevant user manual or other written instructions issued by Nanosonics (including where used in temperature or other external conditions exceeding those set out in the product specification, or serviced by persons other than Nanosonics' approved service personnel);
- b. where modifications have been made to the trophon EPR, other than by Nanosonics or its authorised service providers;
- c. where unauthorised consumables, accessories or other chemicals or items have been used with the trophon EPR;
- d. where the trophon EPR is used in conjunction with other equipment or products (other than multiple use ultrasound probes as described in the user manual), without Nanosonics' prior written consent;
- e. where the trophon EPR has been damaged due to external or environmental causes of any kind (including factors such as voltage fluctuations, excess voltage or power failure);
- f. where the trophon EPR has been damaged as a direct or indirect result of any malicious or negligent act or omission by any person (other than Nanosonics or its authorised service providers);
- g. where the defect does not materially affect the function of the trophon EPR (for example scratches or marks on the external surface of the trophon EPR); or
- h. where the serial number or product label has been removed, changed, deleted or made unrecognizable, or if the number or label is no longer clearly distinguishable for other reasons beyond Nanosonics' control and therefore it is not possible to conclusively identify the product.

This warranty applies to the trophon EPR and defective parts only; the warranty does not cover the replacement of used disinfectant cartridges or of parts which need periodic replacement during the life of the product as a result of the ordinary use made of them.

How to make a claim

Please contact your trophon customer service representative with any queries regarding this warranty or post warranty repairs. If you wish to make a warranty claim, please contact your trophon EPR customer service representative, or contact Nanosonics at [+61 2 8063 1603] or [support@nanosonics.com.au]. You can also send a warranty claim in writing (with proof of purchase) to Nanosonics at: 14 Mars Road, Lane Cove, NSW 2066 Australia.

Nanosonics will make arrangements for and bear the cost of the collection of your trophon EPR. You will be responsible for uninstalling, reinstalling and recommissioning the trophon EPR, regardless of whether or not it is found to be defective. If Nanosonics finds on examination that the trophon EPR is defective in materials and workmanship and is within the warranty period, then we will repair or replace the defective trophon EPR at our discretion. Nanosonics will bear the cost of return delivery of the repaired trophon EPR, or replacement trophon EPR, to you. If we are unable to repair or replace the trophon EPR for any reason, we will discuss with you an appropriate solution including upgrading you to a newer model or refunding the purchase price.

If Nanosonics finds on examination that the trophon EPR is **not** defective in materials and workmanship, or if you are not entitled to the benefit of this warranty (for example, if any of the above

exclusions apply, or the claim was not made within the warranty period), then Nanosonics may require you to bear the cost of returning the trophon EPR to you, and the costs of any repairs to the trophon EPR, or replacement trophon EPR, requested by you.

You are responsible for backing up all data on the trophon EPR if it is being repaired, and acknowledge that repair of the trophon EPR may result in the loss of user-generated data stored on the trophon EPR.

Goods presented for repair may be replaced by refurbished goods of the same type rather than being repaired. Refurbished parts may be used to repair the goods.

Australia: *The following statement only applies if you are a 'consumer' for the purpose of the Australian Consumer Law at Schedule 2 of the Competition and Consumer Act 2010.*

Our goods come with guarantees that cannot be excluded under the *Australian Consumer Law*. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. The benefits under this product warranty are in addition to other rights and remedies you may have under law in relation to our goods.

New Zealand: *The following statement only applies if you are a 'consumer' for the purpose of New Zealand's Consumer Guarantees Act 1993.*

Our goods come with guarantees that cannot be excluded under the *Consumer Guarantees Act 1993*. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

United States: *The following statement only applies to purchasers of the trophon EPR in the United States.*

This warranty is a limited warranty, is the sole and exclusive warranty applicable to the product(s) described herein and is made in lieu of all other warranties, express or implied, including without limitation warranties of merchantability or fitness for a particular purpose.

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Nanosonics Europe GmbH