



CRITIKON SOFT-CUF

Blood Pressure Cuffs

Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It is economically priced to help support your infection control efforts.



Designed for patient safety, comfort and accuracy

- Limited reuse, yet economical enough for pairing with a single patient throughout their hospital stay to help with your infection control efforts
- Soft, absorbent material is not made with DEHP or natural rubber latex to help enhance patient safety
- Designed with rounded corners for patient comfort
- Designed and clinically validated for performance accuracy to help you make well-informed decisions
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit



Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to connect the cuff to the monitor faster than other dual-tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with GE and non-GE automated BP devices, as well as on a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility and optimization of inventory space

Cuff assortment packs available

- For low-volume cuff utilizers, three types of assortment packs are available for easier ordering and more optimal inventory management



| SOFT-CUF | | | ISO 80369-5 compliant | |
|------------------|------------|--------------------|--|---|
| | | |  |  |
| Size | Color | Limb Circumference | 2-Tube DINACLICK Connector | 1-Tube Bayonet Connector |
| Infant | Orange | 8-13 cm | SFT-P1-2A-5S | 2360 |
| Child | Green | 12-19 cm | SFT-P2-2A-5S* | 2361 |
| Child Long | Green | 12-19 cm | SFT-P2-2A-L-5S | |
| Small Adult | Light Blue | 17-25 cm | SFT-A1-2A-5S* | 2362* |
| Small Adult Long | Light Blue | 17-25 cm | SFT-A1-2A-L-5S | |
| Adult | Navy | 23-33 cm | SFT-A2-2A-5S* | 2363* |
| Adult Long | Navy | 23-33 cm | SFT-A2-2A-L-5S* | 2140* |
| Large Adult | Rose | 31-40 cm | SFT-A3-2A-5S* | 2364* |
| Large Adult Long | Rose | 31-40 cm | SFT-A3-2A-L-5S | 2141 |
| Thigh | Brown | 38-50 cm | SFT-T1-2A-5S | 2365 |

Note: 20 cuffs/box

* Has a respective part number with an -INT suffix, denoting that the Instructions For Use (IFU) is available in multiple languages.

DINACLICK Assortment Packs

| Connector Type | Part Number (ISO 80369-5 compliant) | Description |
|-----------------------------|-------------------------------------|--|
| 2-Tube DINACLICK Connectors | 2059303-501 | 5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box) |
| | 2059303-502 | 6 Infant, 10 Child, 4 Small Adult (20 cuffs/box) |
| | 2059303-503 | 3 Small Adult, 7 Adult, 3 Adult Long, 3 Large Adult, 4 RADIAL-CUF (forearm) (20 cuffs/box) |
| | 2059303-504 | 2 Infant, 2 Child, 2 Small Adult, 6 Adult, 2 Adult Long, 2 Large Adult, 2 RADIAL-CUF (forearm), 2 Thigh (20 cuffs/box) |
| | 2059303-506 | 1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF (4 cuffs/pack) |

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at gehealthcare.co.uk/products/patient-monitoring/non-invasive-blood-pressure-cuffs



© 2020 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, DURA-CUF and DINACLICK are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00919XE 8/20