



GE Reusable Bulb Electrode and Limb Clamps

Quality signal transmission.

GE Healthcare reusable bulb electrodes and limb clamps with silver-silver chloride (Ag/AgCl) plating are designed to ensure positive contact with the skin and provide maximum signal quality.

Features

GE Healthcare reusable bulb electrodes and limb clamps are compatible with 4.0 mm banana lead termination leadwires.

Product compatibility

MAC™i, MAC 600, MAC 800, MAC 2000, MAC Link¹ ECG System

Our equipment. Our accessories. Simply compatible.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step

PRODUCT SPECIFICATIONS

Materials

Does not contain natural rubber latex, DEHP, BPA or soft PVC. ROHS (restriction of hazardous substances) and bio-compatibility compliant.

Storage conditions

Recommended temperature -40° to 70°C

Recommended relative humidity 10% to 95%

Cleaning solutions

Ethanol 96%

Hydrogen peroxide 3%

Sodium Fypochlorite 0.5% NaOCl

Isopropyl Alcohol 70%

Mikrozid® sensitive wipes

Kleenaseptic® germicidal and surface cleanser

Cidex Plus® solution

Sporicidin® disinfectant solution

Compliance

AHA and IEC standards

CE Mark

Imagination at work

PART LIST

Part Number	Description
2104783-001	Reusable Bulb Electrode, 6/set
2104784-001	Clamp Limb AHA Value 4/ set
2104785-001	Clamp Limb IEC Value 4/ set

Not available in the US.

Product may not be available in all countries and regions.
Contact a GE Healthcare Representative for more information.
www.gehealthcare.com/promotional-locations

Distributed by:

GE Medical Systems Information Technologies, Inc.



Manufactured by:

Shenzhen Launch Electrical Co., Ltd.
No.F Building, Zhonggangxing
Industrial Estate, Zhangge Community
Guanlan Street, Longhua New District,
Shenzhen, China

© 2018 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or, obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, Imagination at work and MAC are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. Mikrozid is a trademark of Schülke & Mayr GmbH. Kleenaseptic is a trademark and brand of Metrex Research, LLC. Cidex is a trademark of Johnson & Johnson. Sporicidin is a trademark of Contec, Inc. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2180255 12/18